

(19) World Intellectual Property Organization  
International Bureau



(43) International Publication Date  
24 December 2003 (24.12.2003)

PCT

(10) International Publication Number  
**WO 03/106697 A2**

(51) International Patent Classification<sup>7</sup>:

C12Q

(21) International Application Number: PCT/US03/06715

(22) International Filing Date: 19 March 2003 (19.03.2003)

(25) Filing Language: English

(26) Publication Language: English

(30) Priority Data:  
60/365,812 21 March 2002 (21.03.2002) US

(71) Applicant (for all designated States except US): **GONET, Michael, J.** [US/US]; U.S. Department of Veterans Affairs, Office of General Counsel - PSG IV (024), 810 Vermont Avenue N.W., Washington, DC 20420 (US).

(72) Inventor; and

(75) Inventor/Applicant (for US only): **SCIORTINO, Carmen, V., Jr.** [US/US]; U.S. Department of Veterans Affairs, Office of General Counsel - PSG IV (024), 810 Vermont Avenue N.W., Washington, DC 20420 (US).

(74) Agent: **AGARWAL, Dinesh**; Dinesh Agarwal, P.C., 5350 Shawnee Road, Suite 330, Alexandria, VA 22312 (US).

(81) Designated States (*national*): AE, AG, AL, AM, AT, AU, AZ, BA, BB, BG, BR, BY, BZ, CA, CH, CN, CO, CR, CU, CZ, DE, DK, DM, DZ, EC, EE, ES, FI, GB, GD, GE, GH, GM, HR, HU, ID, IL, IN, IS, JP, KE, KG, KP, KR, KZ, LC, LK, LR, LS, LT, LU, LV, MA, MD, MG, MK, MN, MW, MX, MZ, NO, NZ, OM, PH, PL, PT, RO, RU, SC, SD, SE, SG, SK, SL, TJ, TM, TN, TR, TT, TZ, UA, UG, US, UZ, VC, VN, YU, ZA, ZM, ZW.

(84) Designated States (*regional*): ARIPO patent (GH, GM, KE, LS, MW, MZ, SD, SL, SZ, TZ, UG, ZM, ZW), Eurasian patent (AM, AZ, BY, KG, KZ, MD, RU, TJ, TM), European patent (AT, BE, BG, CH, CY, CZ, DE, DK, EE, ES, FI, FR, GB, GR, HU, IE, IT, LU, MC, NL, PT, RO, SE, SI, SK, TR), OAPI patent (BF, BJ, CF, CG, CI, CM, GA, GN, GQ, GW, ML, MR, NE, SN, TD, TG).

**Declaration under Rule 4.17:**

— of inventorship (Rule 4.17(iv)) for US only

**Published:**

— without international search report and to be republished upon receipt of that report

For two-letter codes and other abbreviations, refer to the "Guidance Notes on Codes and Abbreviations" appearing at the beginning of each regular issue of the PCT Gazette.

(54) Title: METHOD AND KIT FOR IDENTIFYING *PSEUDOMONAS AERUGINOSA*

(57) Abstract: A method of identifying a bacteria, such as *Pseudomonas aeruginosa*, in a sample, includes providing a sample suspect of comprising a bacteria to be identified, exposing the sample to an antibody specific for a lipoprotein of the bacteria and an agglutination reagent, allowing the sample to react with the antibody and the agglutination reagent, wherein the presence of the bacteria is indicated if an agglutination occurs. A kit for testing the presence of a bacteria, such as *Pseudomonas aeruginosa*, includes an agglutination reagent and an antibody specific for a lipoprotein of the bacteria.

WO 03/106697 A2

Attorney's Docket No.: PCT 1335/03 (VA)

METHOD AND KIT FOR IDENTIFYING *PSEUDOMONAS AERUGINOSA*

CARMEN V. SCIORTINO, JR., PH.D.

5

CROSS-REFERENCE TO RELATED APPLICATIONS

[0001] The present application claims priority on prior U.S.  
Provisional Application S.N. 60/365,812, filed March 21, 2002, and which is  
10 incorporated herein in its entirety by reference.

BACKGROUND OF THE INVENTION

[0002] The present invention is generally directed to identifying a  
15 bacteria in a sample, and more particularly to a method and kit for identifying  
or assaying *Pseudomonas aeruginosa*.

[0003] *Pseudomonas aeruginosa* is an oxidase-positive, gram-negative, rod-shaped organism that is found ubiquitously in the environment. It is versatile in its habitat and can grow in soil, water, and on plant and animal tissue. It is an opportunistic organism and one of the most problematic nosocomial pathogens capable of causing disease in susceptible individuals such as people who have cystic fibrosis, cancer, burns, or some immune system deficiency (Jaffe, Lane, and Bates 2001). Case fatality can be as high as 50 percent due to a combination of weakened host defenses, bacterial resistance to antibiotics, and the production of extra cellular-bacterial enzymes and toxins (Iglewski 2002). *P. aeruginosa* often colonizes hospital food, sinks, taps, mops, and respiratory equipment. The infection is spread from patient to patient via contact with fomites or by ingestion of contaminated food and water (Iglewski 2002).

15 [0004] *P. aeruginosa* is clinically indistinguishable from other gram-negative bacteria that also cause these sorts of infections but that have a lower morbidity and mortality rate. Therefore, early and accurate diagnosis is important. This is particularly important as *P. aeruginosa* is well known for being resistant to a wide spectrum of antibiotics.

20

[0005] Healthy individuals who come into contact with people with *P. aeruginosa* infections are not at risk of developing the infection themselves (Pseudomonas Genome Project 2002). In fact, *P. aeruginosa* is a

resident of the intestinal tract in about 10 percent of healthy individuals, and is found sporadically in moist areas of the human skin and in the saliva (Chamberlain 2002).

5           **[0006]**           The major clinical features used in diagnosis *in situ* are pus formation, pyocyanin formation in about 90 percent of cases, and fluorescein formation which can be viewed in the dark with a Wood's UV light for fluorescence. There are 13 antigenic groups of *P. aeruginosa*, which in the future, may be treated differentially with immunotherapy (Chamberlain 2002).  
10   This bacterium can be distinguished from other pseudomonad strains by growth at 42°C (Washington State University 2002).

**[0007]**           *P. aeruginosa* is also the most important human pathogen in the genus *Pseudomonas* (Kiska and Gilligan, 1999). It can also  
15   cause superficial skin infections, osteomyelitis, sepsis, folliculitis, swimmer's ear, endocarditis in intravenous-drug users, urinary tract infections, otitis media, and contact lens associated corneal ulcer (Dini et al., 2000; Pollack, 2000). As noted above, *P. aeruginosa* is also a major cause of nosocomial infection (Emori and Gaynes, 1993). In most clinical microbiology  
20   laboratories, *P. aeruginosa* is one of the top three clinical isolates, therefore its frequency of isolation is remarkable. Since *P. aeruginosa* is also highly antibiotic resistant, specific classes of antibiotics are recommended for their anti-pseudomonas activity. Therefore, a 24-48 hour delay in organism

identification may delay appropriate therapy and adversely affect patient morbidity.

[0008] Few individuals outside of the clinical microbiology laboratory setting may realize the difficulty that automated technology has with the identification of *Pseudomonas aeruginosa*. Although automated instruments are advertised as rapid, the best instruments require at least 18 hours to identify *P. aeruginosa*. In contrast, the same instrument can identify *E. coli* in 4 hours. Since a 24 hour pre-growth of the organism is needed, the total average turn-around-time (TAT) for the identification of *P. aeruginosa* is 46-48 hours. Before 1980, the same was true for the identification of *Staphylococcus aureus*, until a rapid latex agglutination test was developed. Now most microbiology laboratories can identify *S. aureus* in 14-18 hours using one of the many rapid tests available. Unfortunately, there has not been forthcoming a rapid test for *P. aeruginosa* until now. The majority (80%) of bacteriological isolates in a microbiology laboratory are one or more of *Staph aureus*, *P. aeruginosa*, and *Escherichia coli*. Therefore, there is a great need in the routine microbiology laboratory for a cost effective rapid test for *P. aeruginosa*.

20

[0009] In 1988, the inventor of the present invention contributed to the development of a monoclonal antibody (MAB PS2) which recognized a 8 KDa outer membrane protein which was later identified as lipoprotein I (LP

l) (Sciortino 1993). LP I was detected by dot blot analysis in the outer membrane of 99.3% of *P. aeruginosa* strains (n=300). Studies were performed comparing the predictive value of MAB PS2 to the biochemical characterization of *Pseudomonas aeruginosa* by two automated instruments, the Vitek Automated Microbiology System (AMS), Sensititre, and a semi-automated instrument, the Dynatech 2000 (Sciortino et al. 1987). The evaluation showed an 89.7% agreement on organism identification between instruments. Discrepant results between MAB PS2 and the Dynatech 2000, the Sensititre, or the AMS were checked by conventional tubed media identification schemes or the VAMC reference laboratory. In the screening of clinical isolates reported as *P. aeruginosa* by the Dynatech 2000, MAB PS2 negatively reacted with two organisms which were later determined by conventional tests not to be *P. aeruginosa* since they did not grow on centrimide or at 42°C.

15

[0010] In screening several stock cultures of gram-negative organisms, two strains which were previously unidentified and listed only by a catalog number gave positive reactions with PS2. These strains grew on centrimide and at 42°C which is indicative of *P. aeruginosa*. When tested by the AMS, one strain was identified as *P. aeruginosa* with 99% confidence while the other strain only showed a 14% probability of *being P. aeruginosa*. Because this organism grew on centrimide and at 42°C, it was correctly

identified by MAB PS2 as *P. aeruginosa*. This indicates that the MAB PS2 can be more sensitive in the recognition of this organism than the AMS.

**[0011]** Recent studies have involved development of test  
5 reagents and definition of specific reactivities of the molecules. For instance, the inventor discovered that *P. aeruginosa* (LP I) expression is time dependent. Twenty-four hour cultures express very little LP I in the outer membrane whereas 48 h cultures have large amounts of LP I associated with the outer membrane. It then became necessary to extract LP I from the inner  
10 membrane for maximum antibody reactivity at 24 h. Another problematic technical difficulty in the past studies involved adjusting trace metal concentrations. Too much magnesium caused auto-agglutination of bacteria with the reagents although  $Mg^{++}$  was required to release LP I from the bacterial membrane complex. Therefore, recent studies have focused on  
15 optimization of the co-agglutination test and less on the sensitivity and specificity of the test.

**[0012]** Diagnosis of *P. aeruginosa* therefore, depends on its isolation and laboratory identification. It can be cultured on most general-  
20 purpose media and is commonly isolated on blood agar plates or eosin-methylthionine blue agar. It is identified on the basis of its gram morphology, inability to ferment lactose, a positive oxidase reaction, its fruity odor, and its

ability to grow at 42° C. It fluoresces under ultraviolet radiation and this is useful in suggesting its presence in wounds (Iglewski 2002).

[0013] The current automated technology for identification of *P. aeruginosa* is Vitek Automated Microbiology System. This method costs approximately \$7 per sample (Jaffe, Lane, and Bates, 2001).

[0014] Biotec Laboratories Limited, a U.K. biotechnology company, is focusing on a novel technology termed phage amplification that is based on the presence of a specific bacteriophage for identifying the presence of the target bacteria within a 4-hour period (Stewart et al. 1998). Biotec's phage amplification is a platform technology with a test for *P. aeruginosa* in the research and development phase (per personal communication). The company claims this technology will have the advantages of speed, accuracy, simplicity, and low cost (Biotec Laboratories Limited 2002).

[0015] Jaffe, Lane, and Bates (2001) reportedly have developed a rapid DNA extraction and polymerase chain reaction (PCR) method for identifying *P. aeruginosa* directly from clinical samples including blood, sputum, and urine. Two versions of this test have been developed: standard PCR and rapid-cycle PCR, with the former having a turnaround time of 4 hours and the latter of just 1 hour. The results obtained using this technique



are 100 percent in concordance with those achieved using the conventional culture methods. They estimate the direct supply cost of the DNA extraction and PCR to be \$2.22 per sample. It is important to note, however, that the equipment costs involved in real-time PCR technology are considerable and have almost certainly not been factored into this calculation. This work has been accomplished in collaboration with Commonwealth Biotechnologies, Inc. (Virginia), a contract research organization. The developers do not appear to be selling this kit, nor does Commonwealth Biotechnologies provide any further information about it on the company's web site (Commonwealth Biotechnologies, Inc. 2002).

**[0016]** Therefore, there has not yet been a rapid test for *P. aeruginosa* that would allow its identification on the first day of culture. Automated microbiology instruments still require 18-24 h to identify *Pseudomonas*. Roche Inc, has recently introduced a rapid light-cycler test for rapid identification of *P. aeruginosa*, but equipment cost and test complexity may limit its use in the routine clinical laboratory.

**[0017]** In view of the above, there is a need in the industry for a rapid, accurate, reliable and inexpensive test for the identification of *Pseudomonas aeruginosa*.

## OBJECTS AND SUMMARY OF THE INVENTION

[0018] The principal object of the present invention is to provide a method and kit for identifying *Pseudomonas aeruginosa* which overcomes  
5 the drawbacks associated with conventional techniques and materials.

[0019] An object of the present invention is to provide a method and kit for identifying *Pseudomonas aeruginosa* which is fast, reliable, accurate, sensitive, specific, and inexpensive.

10

[0020] Another object of the present invention is to provide a method and kit for identifying *Pseudomonas aeruginosa* which significantly reduces the test turnaround time to about 18-20 hours, and is a single assay as opposed to a battery of tests required by the conventional methods. In  
15 particular, the instant invention reduces the laboratory analysis turnaround time by about 24 hours, providing diagnosis the day after the sample is received and consequently allowing antibiotic or other suitable treatment to begin a day earlier. This is significant given the morbidity associated with this infection. The test has a sensitivity of about 99.3% and a specificity of about  
20 95%.

[0021] Yet another object of the present invention is to provide a method and kit for identifying *Pseudomonas aeruginosa* which uses an

antibody that is specific for a lipoprotein (LP1) on the surface of *Pseudomonas aeruginosa*, and an agglutination reagent to indicate a positive result.

5           **[0022]**       A further object of the present invention is to provide a method and kit for identifying *Pseudomonas aeruginosa* which costs significantly less than the conventional tests.

**[0023]**       Yet a further object of the present invention is to provide  
10   a method and kit for identifying *Pseudomonas aeruginosa* which is more likely to be accepted by the clinical diagnostics community because it is an antibody-based diagnostic test for bacterial infections than the conventional tests based on other formats, such as DNA diagnostics.

15           **[0024]**       Still a further object of the present invention is to provide a method and kit for identifying *Pseudomonas aeruginosa* which would fill a market need for a rapid and reliable test for a problematic pathogen that causes significant morbidity and mortality.

20           **[0025]**       An additional object of the present invention is to provide a method and kit for identifying *Pseudomonas aeruginosa* which rapidly identifies *Pseudomonas aeruginosa* and distinguishes it from other microorganisms. Rapid testing is important because *Pseudomonas*

*aeruginosa* is highly resistant to most antibiotics and a faster identification translates into faster suitable treatment.

[0026] In summary, the main object of the present invention is to  
5 provide a method and kit for identifying *Pseudomonas aeruginosa* which is fast, reliable, accurate, sensitive, specific and inexpensive. The test of the invention would result in faster treatment of a bacterial infection due to the rapid identification of *Pseudomonas aeruginosa* by the technique of the present invention.

10

[0027] In accordance with the present invention, a method of identifying a bacteria in a sample, includes providing a sample suspect of comprising a bacteria to be identified, exposing the sample to an antibody specific for a lipoprotein of the bacteria and an agglutination reagent, allowing  
15 the sample to react with the antibody and the agglutination reagent, wherein the presence of the bacteria is indicated if an agglutination occurs.

[0028] In accordance with the present invention, a kit for testing the presence of a bacteria, such as *Pseudomonas aeruginosa*, in a sample,  
20 includes an agglutination reagent and an antibody specific for a lipoprotein of the bacteria.

## DETAILED DESCRIPTION OF THE INVENTION

[0029] The present invention is directed to a method or test (named "Pseudostat II") for identifying *P. aeruginosa* from first day cultures with only minimum alternative testing and equipment that is common to most clinical laboratories. The basis for the test is a monoclonal antibody (Mab) that targets lipoprotein I which is unique to *P. aeruginosa* (De Vos et al., 1993; De Vos et al., 1998; Mutharia et al., 1982; Saint-Onge et al., 1992).

10 [0030] In carrying out the present invention, bacteria were cultivated using standard clinical microbiology procedures, with incubation at 37°C, 5% CO<sub>2</sub>, on trypticase-soy agar with 5% sheep blood. Stock quality control strains and some of our immuno-type strains were stored at -70°C as suspensions in BHI-broth, 20% glycerol. The clinical isolates used in this study were freshly obtained in-house from patients with gram-negative bacterial infections.

[0031] Lipoprotein I specific Mab PS2, was prepared by growth of hybridomas in tissue culture flasks using serum and protein free medium (Sigma, St. Louis, MO). Mab was purified by ammonium sulfate precipitation followed by dialysis against 0.1 M Trizma. *Staphylococcus aureus* reagent was prepared by formalin treatment of 24 h cultures of *S. aureus* (ATCC # 12598, Cowan serotype I), followed by numerous centrifugations and

washings in 0.1 M Trizma buffer, pH 7.2, and staining with methylene blue. Extraction and reactivity buffers were prepared using commercially available reagents.

5           **[0032]**       The Pseudostat II test kit included four reagents. Namely, *S. aureus*, "reagent A"; buffered Mab PS2, "reagent B"; negative control reagent, "reagent C", and "extraction buffer". An oxidase spot-test was performed on all 18 and 24 h old isolates. Only oxidase-positive isolates were tested with the Pseudostat II test. All isolates were previously identified using  
10 the Vitek gram-negative identification card. Low probability identification calls (<85%) by the Vitek were additionally tested using the API 20 NE test system (bioMérieux Vitek, Inc. Hazelwood, MO, USA.). Additional tests used to verify isolate identity were: growth on Pseudosel agar (Becton Dickinson Microbiology Systems, Cockeysville, Maryland, USA), growth at 42°C, and  
15 fluorescence under ultraviolet light when grown on Mueller Hinton agar (Kiska and Gilligan, 1999).

**[0033]**       The test procedure was performed by removing one or two colonies from agar plates using the cotton tip end of a sterile cotton-tip  
20 applicator stick. The swab containing bacteria was immersed into 0.5 ml of extraction buffer (in a 1 ml conical microfuge tube), then twisted and compressed against the side of the tube. The tube was capped and vortexed for 10 sec and then placed in a 37°C water bath for 5 min. White Time tape,

3/4" width by 1.5" length, was stuck to an index card for use as the agglutination card. At the end of the extraction, samples were again vortexed for 10 sec. Ten  $\mu$ l of sample were removed with a disposable pipet tip and mixed on the agglutination card with 1 drop of reagent A, and 1 drop of reagent B, then spread in a 2-cm<sup>2</sup>-diameter oval. The index card with reagents was placed on a lab rotator for 5-8 min. A positive test was recognized as a blue agglutination against a white background. Three controls were run each date of testing. The positive control was *P. aeruginosa* ATCC #27853. The negative control was *Bordetella bronchiseptica* ATCC #10580. The reagent control consisted of the Mab suspension buffer without the Mab. A negative test was recognized by the absence of agglutination or agglutination equivalent to that of the reagent control, scored -, or  $\pm$ , respectfully. Positives were scored as 1+ for granular agglutination, 2+ for flocculate agglutination, 3+ if a string formed, and 4+ if a button formed. Strong 4+ agglutination occurred within 3-5 min whereas weaker reactions took 4-8 min.

[0034] A total of 232 *P. aeruginosa* isolates and 36 other oxidase-positive, gram-negative clinical isolates were tested. Included in the *P. aeruginosa* isolates were the International Antigenic Type (IAT) strains and 6 mucoid clinical isolates. 75% of the *P. aeruginosa* isolates including all of the mucoid and IAT strains gave 4+ reactions (Table 1 below). The sensitivity of the test was 96% and the specificity of the test was 91.9% (Table 2 below).

Only nine isolates were falsely negative, however one of these, when retested gave a 3+ reaction. In this study, only the first reaction was used with the consideration that if this were being performed in a clinical laboratory, only one test would be performed. The three false-positives, *Achromobacter* 5 *xylosoxidans*, *Burkholderia cepacia*, and *Pseudomonas fluorescens*, were repeated but consistently gave false-positive reactions. All were re-tested using the API 20 NE system and all gave >95% probabilities for identification. The *P. fluorescens* isolate was a mucoid environmental isolate that grew at room temperature but not at 37°C. Although it was an exception to study 10 criteria, its inclusion showed that some false-positives occurred if limitations of the test were exceeded.

[0035] The reactivity of Mab PS2, has been previously investigated for its reactivity against 470 bacterial isolates using the bio-dot 15 technique (Sciortino, 1993). In that study, Mab PS2 recognized 298/300 *P. aeruginosa* isolates that included all IAT strains and the Fisher-Devlin immuno-type strains. It did not react with 118/123 other gram-negative bacteria. The outliers in that study to which, Mab PS2 reacted were 3/17 *Enterobacter aerogenes*, 1/11 *Escherichia coli*, and 1/2 *Proteus vulgaris* 20 isolates. For this reason, the use of this test is preferably limited to only oxidase-positive bacteria. Another preference for the test of the invention is that cultures be incubated for at least 18 h. This is because young cultures may not express sufficient amounts of lipoprotein I on the cell surface. For



example, isolates that gave 1+ to 3+ reactions when incubated an additional 24 h, all gave 4+ reactions. Some bacteria that were tested showed weak reactivity with the reagent control. This was attributed to the non-specific binding of *S. aureus* with some component of the digested bacterial mixture.

- 5 The reagent-control reactivity was subtracted from the isolate-test reactivity for corrected interpretation of reactions.

[0036] The following are the details for the method/test of the present invention.

10

#### Reagent Preparation

- [0037] Reagent A (Staph-A co-agglutination reagent). The Cowan strain of *Staphylococcus aureus* was purchased from the American Type Culture Collection. The reagent was prepared by formalin treatment of 24 h cultures of *Staphylococcus aureus* followed by numerous washings in buffer and staining with methylene blue.
- 15

20

[0038] Staph A Reagent included the following:

- Stock reagent of Formalin killed, methylene-blue dyed

*Staphylococcus aureus* in 0.1 M HEPES pH 7.6

-500 mM MES buffer, pH 6.0- Suspension buffer

5 -10% FSG (Fish Skin Gelatin) – non-specific blocking  
reagent

-Sodium azide - preservative

[0039] Reagent B (Monoclonal antibody). Monoclonal antibody

10 PS2 was purchased commercially. (IMMR, Charleston, S. Carolina. Contact  
Person: Amy Reid 118 St. Michaels Place, Monks Corner, S.C. 29461)

[0040] PS2 Reagent B included the following:

-Mab PS2: 3 mg/ml of 0.01 M Trizma buffer, pH 7.6 –

15 Stock antibody

-500 mM MES buffer, pH 6.0 – Stock suspension buffer

-0.2 M EDTA [(ethylenedinitrillo)- tetraacetic acid] – Stock  
chelating agent

-10% FSG (Fish Skin Gelatin) in water - Stock non-  
20 specific blocking reagent Lysozyme , 10 mg/ml

dissolved in 0.01 M glycine/Trizma/NaCl pH 2.5 -

Stock digestion enzyme

-Sodium azide - preservative

**[0041]** Reagent C (Negative control reagent) included the following:

- 500 mM MES buffer, pH 6.0 – Stock suspension buffer
- 0.2 M EDTA- Stock chelating agent
- 10% FSG (Fish Skin Gelatin) in water - Stock non-specific blocking reagent
- Sodium azide - preservative

**[0042]** Reagent D (Extraction buffer). Extraction and reactivity buffers are prepared in-house using commercially available reagents.

**[0043]** Extraction Buffer included the following:

- 500 mM MES pH 6.0 - Stock suspension buffer
- 1% Triton X 100 in water – Stock suspension detergent
- Lysozyme Stock: 10 mg/ml dissolved in 0.01 M glycine/Trizma/NaCl pH 2.5
- MgCl<sub>2</sub>, 0.1M - Stock
- Sodium azide- preservative

User Kit

**[0044]** The kit of the invention includes the above-identified three reagents and the negative control reagent. The negative control reagent C, is

the same as Reagent B, but lacks the monoclonal antibody. Each reagent is labeled, A,B,C, and Extraction buffer. Preferably, kits contain 100 tests per kit. The kit also contains micro-centrifuge tubes, cotton swabs, tube rack, and time tape with index cards.

5

## TEST PROCEDURE

**[0045]** Preferably, bacterial cultures must be 18 hours old. An oxidase spot-test is performed on suspected colonies. Only oxidase-positive isolates are tested further. The oxidase test is a standard microbiological test performed in all Microbiology laboratories. It requires about two minutes of test time.

### Step One

15 **[0046]** Oxidase-positive colonies are removed from agar plates using a sterile cotton-tip applicator stick. A generous portion of culture is used, i.e., 2-5 colonies. The cotton tip and not the wood end of the applicator stick is used.

### 20 Step Two

**[0047]** The cotton swab containing bacteria is immersed into 0.5 ml of extraction buffer (in a microfuge tube).

**Step Three**

**[0048]** Twist the swab in the tube and wring the swab against the sides of the tube.

**5 Step Four**

**[0049]** Cap the tube and place @ 37°C for 5-10 minutes. (A water bath or heat block is preferable.)

**Step Five**

**10 [0050]** Place a white piece of Time tape, preferably 3/4" wide and 1.5" in length, on an index card.

**Step Six**

**[0051]** At the end of the extraction period, with the tube capped, briefly vortex the sample. Then mix reagents on the piece of Time tape as follows:

-1 drop of Staph A (Reagent A)

-1 drop of Antibody reagent (Reagent B)

-Add 10 µl of extracted bacterial suspension from the Extraction buffer mix.

**20**

### Step Seven

**[0052]** Mix the reagents with the pipet tip, forming an oval circle about the size of a nickel.

### 5 Step Eight

**[0053]** Place the index card with reagents on a lab rotator and let rotate for about 5 minutes.

## INTERPRETATION OF THE TEST RESULTS

10

**[0054]** A positive test is recognized by its agglutination reaction. It is graded, from zero to 4+ reaction. The agglutination reaction is a blue flocculate precipitate against a white background. Some auto-agglutination of reagents may occur for some strains of bacteria. The agglutination seen in  
15 the test sample must be greater than the negative control to establish a positive test result. (See the Negative Control procedure below.)

**[0055]** A negative test is recognized by the absence of agglutination.

20

Negative Control (Performed simultaneously with the test organism above.)

**[0056]** For each organism tested, a negative reagent control is run as follows:

### Step One

**[0057]** Place a white piece of Time tape, preferably 3/4" wide and 1.5" in length, on an index card next to the piece of tape used for the  
5 primary test above.

### Step Two

**[0058]** At the end of the extraction period, briefly vortex the reaction tube. Then mix reagents on the piece of Time tape as follows:  
10 -1 drop of Staph A (Reagent A)  
-1 drop of reagent control (Reagent C)

**[0059]** Mix the reagents with the pipet tip, forming an oval circle about the size of a nickel.  
15

**[0060]** Place the index card with reagents on a lab rotator and let it rotate for 5 minutes.

**[0061]** A negative test is recognized by the absence of flocculate  
20 agglutination.

### Quality Control

[0062] Quality control is performed with a strain of *Pseudomonas aeruginosa* (ATCC #27853) and with *Bordetella bronchioseptica* ATCC strain, *P. aeruginosa* (ATCC #27853) gives a positive  
5 reaction, whereas *B bronchioseptica* gives a negative reaction.

[0063] From the above, one can readily observe that the Pseudostat II co-agglutination test of the invention would be particularly important to clinical microbiology laboratories for the first day identification of  
10 *P. aeruginosa*. When used in conjunction with colonial morphology, odor, pigmentation, and oxidase-reactivity, it would provide an accurate and rapid means to identify *P. aeruginosa* and to alert physicians early that *P. aeruginosa* specific antibiotics should be considered for therapy.

15

20



TABLE 1  
Reactivity of the Pseudostat II Test for Oxidase-Positive,  
Gram-Negative Isolates

Organism	No. of Isolates	Reactivity	Interpretation
<u>Achromobacter xylosoxidans</u>	3	-	Negative
<u>subsp. denitrificans</u>			
<u>Achromobacter xylosoxidans</u>	1	4+	Positive
<u>subsp. xylosoxidans</u>			
<u>Achromobacter xylosoxidans</u>	2	-	Negative
<u>subsp. xylosoxidans</u>			
<u>Alcaligenes odorans</u>	1	-	Negative
<u>Bordetella bronchiseptica</u>	2	-	Negative
<u>Burkholderia cepacia</u>	2	-	Negative
<u>Burkholderia cepacia</u>	1	2+	Positive
<u>*Pseudomonas oryzihabitans</u>	2	-	Negative
<u>Moraxella catarrhalis</u>	8	-	Negative
<u>Moraxella lacunata</u>	1	-	Negative
<u>Pseudomonas aeruginosa</u>	2	-	Negative
<u>Pseudomonas aeruginosa</u>	7	±	Negative
<u>Pseudomonas aeruginosa</u>	5	1+	Positive
<u>Pseudomonas aeruginosa</u>	20	2+	Positive
<u>Pseudomonas aeruginosa</u>	23	3+	Positive
<u>Pseudomonas aeruginosa</u>	175	4+	Positive
<u>Pseudomonas fluorescens</u>	1	-	Negative
<u>Pseudomonas fluorescens</u>	1	2+	Positive
<u>Pseudomonas stutzeri</u>	1	-	Negative
<u>Pseudomonas vesicularis</u>	1	-	Negative
<u>Pseudomonas putida</u>	1	-	Negative
<u>Ralstonia pickettii</u>	3	-	Negative
<u>Ralstonia pickettii</u>	1	±	Negative
<u>Shewanella putrefaciens</u>	1	-	Negative
<u>Shewanella putrefaciens</u>	1	±	Negative
<u>Sphingomonas paucimobilis</u>	3	-	Negative

5 \*Weakly oxidase positive isolates

-, negative with no agglutination; ±, negative but equal to the reagent negative control; 1+, positive granular agglutination; 2+ positive flocculate agglutination; 3+, positive string formation with agglutination; 4+, positive agglutination with button formation

10

TABLE 2  
Sensitivity and Specificity of the Pseudostat II Test

Group	No. of Isolates		Sensitivity of test (%)	Specificity of test (%)
	Positive	Negative		
<i>Pseudomonas aeruginosa</i>	223	9	96	91.9
Other gram-negatives	3	34		

5

10

15

20

[0064] While this invention has been described as having preferred sequences, ranges, steps, materials, or designs, it is understood that it includes further modifications, variations, uses and/or adaptations thereof following in general the principle of the invention, and including such  
5 departures from the present disclosure as those come within the known or customary practice in the art to which the invention pertains, and as may be applied to the central features hereinbefore set forth, and fall within the scope of the invention and of the limits of the appended claims.

10 [0065] The following references, to the extent that they provide exemplary procedural or other details supplementary to those set forth herein, are specifically incorporated herein by reference.

#### REFERENCES

15

Biotec Laboratories Limited. 2002. Biotec In Development. Web site:

<http://www.biotec.co.uk/development.html>. (April 29, 2002.)

BIVDA. 2000. British In Vitro Diagnostics Association. IVD - Healthcare Survey.

Web site: <http://www.bivda.co.uk/laboratory/index.html>. (April 16, 2002.)

20 Chamberlain, N. 2002. Pseudomonas infections. Kirksville College of

Osteopathic Medicine. Web site:

<http://www.kcom.edu/faculty/chamberlain/Website/tritzid/pseudinf.htm>

(April 29, 2002.)

- Commonwealth Biotechnologies, Inc. 2002. Web site: <http://www.cbi-biotech.com>. (April 29, 2002.)
- De Vos, D., Bouton, C., Sarniguet, A., De Vos, P., Vauterin, M., and P. Cornelis. 1998. Sequence diversity of the *oprI* gene, coding for major outer membrane lipoprotein I, among rRNA group I Pseudomonads. J. Bacteriol. 180: 6551-6556.
- De Vos, D., Lim, A. Jr., De Vos, P., Sarniguet, A., Kersters, K., and P. Cornelis. 1993. Detection of the outer membrane lipoprotein I and its gene in fluorescent and non-fluorescent pseudomonads: implications for taxonomy and diagnosis. J. Gen. Microbiol. 139:2215-2223.
- Dini, L.A., Cockinos, C., Frea, J.A., Niszl, I.A., and M.B. Markus. 2000. Unusual case of *Acanthamoeba polyphaga* and *Pseudomonas aeruginosa* keratitis in a contact lens wearer from Gauteng, South Africa. J. Clin. Microbiol. 38:826-829.
- Emori, T.G. and R.P. Gaynes. 1993. An overview of nosocomial infections, including the role of the microbiology laboratory. Clin. Microbiol. Rev. 10:428-442.
- Finke, M., Duchene, M., Eckhardt, A., Domdey, H., and B.U. von Specht. 1990. Protection against experimental *Pseudomonas aeruginosa* infection by recombinant *P. aeruginosa* lipoprotein I expressed in *Escherichia coli*. Infect. Immun. 58:2241-2244.

- Iglewski, B. 2002. *Pseudomonas*. UTMB Graduate School of Biomedical Sciences. Web site: <http://gsbs.utmb.edu/microbook/ch027.htm>. (April 29, 2002.)
- 5 Jaffe, R.I., J.D. Lane, and C.W. Bates. 2001. Real-time identification of *Pseudomonas aeruginosa* direct from clinical samples using a rapid extraction method and polymerase chain reaction. *Journal of Clinical Laboratory Analysis* 15:131-137.
- Kiska, D.L., and P.H. Gilligan. 1999. *Pseudomonas*. In: Murray, P.R., Baron, E.J., Tenover, F.C., and R. H. Tenover (Eds.), *Manual of Clinical Microbiology*, 7<sup>th</sup> ed. American Society for Microbiology, 10 Washington, D.C., USA, pp. 517-525.
- Mutharia, L.M, Nicas, T.I., and R.E.W. Hancock. 1982. Outer membrane proteins of *Pseudomonas aeruginosa* serotype strains. *J. Inf. Dis.* 146:770-779.
- 15 Pollack, M. 2000. *Pseudomonas aeruginosa*. In: Mandell, G. L., Bennett, J.E., and R. Dolin (Eds.). *Principles and practice of infectious diseases*, 5<sup>th</sup> ed. Churchill Livingstone, Philadelphia, Pennsylvania, USA, pp. 2310-2335.
- Pseudomonas Genome Project. 2002. *Pseudomonas aeruginosa* Community 20 Annotation Project. Web site: <http://www.cmdr.ubc.ca/bobh/PAAP.htm> (April 29, 2002.)
- Saint-Onge, A., Romeyer, F., Lebel, P., Masson, L. and R. Brousseau. 1992. Specificity of the *Pseudomonas aeruginosa* PAO1 lipoprotein I gene as

a DNA probe and PCR target region within the *Pseudomonadaceae*. J. Gen. Microbiol. 138:733-341.

Sciortino, C.V. 1993. Surveillance of gram-negative bacteria for *Pseudomonas aeruginosa* lipoprotein I using a monoclonal antibody. Hybridoma 12:327-332.

Sciortino, C.V., 1993. Surveillance of Gram-negative Bacteria for *Pseudomonas aeruginosa* Lipoprotein I using a Monoclonal Antibody. Hybridoma , 12 (3):333-342.

Sciortino, C.V., and A.C.Holm. 1987. Monoclonal antibodies to *Pseudomonas aeruginosa* aminoglycoside-resistance factors. 87th Annual Meeting of the American Society of Microbiology. pp.34, B56.

Stewart G.S.A.B., Jassim, S.A.A., Denyer, S.P., Newby, P., Linley, K. and Dhir, V.K. 1998. The specific and sensitive detection of bacterial pathogens within 4 hours using bacteriophage amplification. *Journal of Applied Microbiology* 84, 777-783.

U.S. Food and Drug Administration Center for Devices and Radiological Health. 2002 . An Introduction to Medical Device Regulations (FDA 92-4222)

Washington State University. Biology Bulletin Board. *Pseudomonas*. 2002. Web site: <http://www.sci.wsu.edu/bio/micro310pseudomonas.html> (April 29, 2002.)

## WHAT IS CLAIMED IS:

1. A method of identifying a bacteria in a sample, comprising the steps of:
  - a) providing a sample suspect of comprising a bacteria to be identified;
  - 5 b) exposing the sample to an antibody specific for a lipoprotein of the bacteria and an agglutination reagent;
  - c) allowing the sample to react with the antibody and the agglutination reagent; and
  - 10 d) whereby the presence of the bacteria is indicated if an agglutination occurs.
2. The method of Claim 1, wherein:

the bacteria to be identified comprises an oxidase-positive organism.
3. The method of Claim 2, wherein:

the antibody comprises monoclonal antibody PS2.

4. The method of Claim 3, wherein:  
the agglutination reagent comprises a strain of *Staphylococcus* bacteria.
5. The method of Claim 4, wherein:  
the organism comprises *Pseudomonas aeruginosa*.
6. A method of identifying *Pseudomonas aeruginosa* in a sample, comprising the steps of:
- a) providing a sample suspect of comprising *Pseudomonas aeruginosa*;
  - b) cultivating organism in a suitable medium;
  - c) removing a portion of the cultured organism and exposing to a first reagent;
  - d) exposing the organism obtained in step c) to second and third reagents, the second reagent comprising an agglutination reagent and the third reagent comprising an antibody specific for a lipoprotein of *Pseudomonas aeruginosa*;
  - e) allowing the components in step d) to react; and
  - f) whereby the presence of *Pseudomonas aeruginosa* is indicated if an agglutination occurs.



5

7. The method of Claim 6, wherein:  
the step c) of removing and exposing a portion of the cultured organism to a first reagent comprises extracting from the *Pseudomonas aeruginosa* organism a lipoprotein by an extraction reagent.
8. The method of Claim 7, wherein:  
the lipoprotein comprises lipoprotein 1.
9. The method of Claim 6, wherein:  
the agglutination reagent comprises a strain of *Staphylococcus* bacteria.
10. The method of Claim 9, wherein:  
the *Staphylococcus* strain is stained generally blue; and  
the presence of *Pseudomonas aeruginosa* is indicated by a generally blue agglutination.
11. The method of Claim 9, wherein:  
the antibody comprises monoclonal antibody PS2.

12. The method of Claim 6, wherein:  
the step b) comprises incubating for a period of about 18 to 24 hours.
13. The method of Claim 6, further comprising the step of:  
g) isolating oxidase-positive organism from the medium prior to step c).
14. A kit for testing the presence of *Pseudomonas aeruginosa* in a sample, comprising:  
a) an agglutination reagent; and  
b) an antibody specific for a lipoprotein of *Pseudomonas aeruginosa*.
15. The kit of Claim 14, further comprising:  
a) a reagent for extracting the lipoprotein from *Pseudomonas aeruginosa*.
16. The kit of Claim 15, wherein:  
a) the agglutination reagent comprises a strain of *Staphylococcus* bacteria.

17. The kit of Claim 16, wherein:
  - a) the antibody comprises monoclonal antibody PS2.
  
18. The kit of Claim 15, further comprising:
  - a) a negative control reagent.

(19) World Intellectual Property  
Organization  
International Bureau



(43) International Publication Date  
24 December 2003 (24.12.2003)

PCT

(10) International Publication Number  
**WO 2003/106697 A3**

(51) International Patent Classification<sup>7</sup>: **G01N 33/53**,  
33/567, C07K 16/00

(21) International Application Number:  
PCT/US2003/006715

(22) International Filing Date: 19 March 2003 (19.03.2003)

(25) Filing Language: English

(26) Publication Language: English

(30) Priority Data:  
60/365,812 21 March 2002 (21.03.2002) US

(71) Applicant (for all designated States except US): **US DEPARTMENT OF VETERANS AFFAIRS** [US/US]; Office of General Counsel - PSG IV (024), 810 Vermont Avenue N.W., Washington, DC 20420 (US).

(72) Inventor; and

(75) Inventor/Applicant (for US only): **SCIORTINO, Carmen, V., Jr.** [US/US]; U.S. Department of Veterans Affairs, Office of General Counsel - PSG IV (024), 810 Vermont Avenue N.W., Washington, DC 20420 (US).

(74) Agent: **AGARWAL, Dinesh**; Dinesh Agarwal, P.C., 5350 Shawnee Road, Suite 330, Alexandria, VA 22312 (US).

(81) Designated States (national): AE, AG, AL, AM, AT, AU, AZ, BA, BB, BG, BR, BY, BZ, CA, CH, CN, CO, CR, CU,

CZ, DE, DK, DM, DZ, EC, EE, ES, FI, GB, GD, GE, GH, GM, HR, HU, ID, IL, IN, IS, JP, KE, KG, KP, KR, KZ, LC, LK, LR, LS, LT, LU, LV, MA, MD, MG, MK, MN, MW, MX, MZ, NO, NZ, OM, PH, PL, PT, RO, RU, SC, SD, SE, SG, SK, SL, TJ, TM, TN, TR, TT, TZ, UA, UG, US, UZ, VC, VN, YU, ZA, ZM, ZW.

(84) Designated States (regional): ARIPO patent (GH, GM, KE, LS, MW, MZ, SD, SL, SZ, TZ, UG, ZM, ZW), Eurasian patent (AM, AZ, BY, KG, KZ, MD, RU, TJ, TM), European patent (AT, BE, BG, CH, CY, CZ, DE, DK, EE, ES, FI, FR, GB, GR, HU, IE, IT, LU, MC, NL, PT, RO, SE, SI, SK, TR), OAPI patent (BF, BJ, CF, CG, CI, CM, GA, GN, GQ, GW, ML, MR, NE, SN, TD, TG).

**Declaration under Rule 4.17:**

— of inventorship (Rule 4.17(iv)) for US only

**Published:**

— with international search report  
— before the expiration of the time limit for amending the claims and to be republished in the event of receipt of amendments

**(88) Date of publication of the international search report:**

19 February 2004

For two-letter codes and other abbreviations, refer to the "Guidance Notes on Codes and Abbreviations" appearing at the beginning of each regular issue of the PCT Gazette.

(54) Title: METHOD AND KIT FOR IDENTIFYING *PSEUDOMONAS AERUGINOSA*

(57) Abstract: A method of identifying a bacteria, such as *Pseudomonas aeruginosa*, in a sample, includes providing a sample suspect of comprising a bacteria to be identified, exposing the sample to an antibody specific for a lipoprotein of the bacteria and an agglutination reagent, allowing the sample to react with the antibody and the agglutination reagent, wherein the presence of the bacteria is indicated if an agglutination occurs. A kit for testing the presence of a bacteria, such as *Pseudomonas aeruginosa*, includes an agglutination reagent and an antibody specific for a lipoprotein of the bacteria.

# INTERNATIONAL SEARCH REPORT

International application No.

PCT/US03/06715

<b>A. CLASSIFICATION OF SUBJECT MATTER</b> IPC(7) : G01N 33/53,567; C07K 16/00 US CL : 435/7.1,7.2; 530/387.1,387.5,388.1,388.2 According to International Patent Classification (IPC) or to both national classification and IPC		
<b>B. FIELDS SEARCHED</b> Minimum documentation searched (classification system followed by classification symbols) U.S. : 435/7.1,7.2; 530/387.1,387.5,388.1,388.2 Documentation searched other than minimum documentation to the extent that such documents are included in the fields searched Bergey's manual Electronic data base consulted during the international search (name of data base and, where practicable, search terms used) Dialog, EAST,		
<b>C. DOCUMENTS CONSIDERED TO BE RELEVANT</b>		
Category *	Citation of document, with indication, where appropriate, of the relevant passages	Relevant to claim No.
Y	US 5,869,272 A (BOGART et al) 9 February 1999 (09.02.1999), column 2, lines 25-47.	1-18
X	US 5,716,829 A (ROSOK et al) 10 February 1998 (10.02.1998), column 4, lines 40-41; column 6, lines 22-42.	2, 5, 6, 7, 14-15, 18
T	US 2003/0170613 A1 (STRAUS) 11 September 2003, entire article especially pages 1, paragraph, 14; page 6, paragraph 69, page 5, paragraph 62, and 16-17.	1-2, 5, 6, 9, 10
---		-----
Y		3-4, 7-8, 11
X	US 5,635,596 A (CHAMBON et al) 3 June 1997 (03.06.1997), column 1, line 10-22 and Fig. 1.	3, 11
---		-----
Y		1-2, 4-10, 12-18
X	SCIORTINO, JR., C.V. Surveillance of gram-negative bacteria for Psuedomonas aeruginosa lipoprotein I using a monoclonal antibody. Hybridoma, 1993, Vol. 12, No. 3, pages 327-332, entire article.	1-3, 8, 11-13, 5-8, 11-15, 17-18
---		-----
Y		4, 9-10, 16
<input checked="" type="checkbox"/> Further documents are listed in the continuation of Box C. <input type="checkbox"/> See patent family annex.		
* Special categories of cited documents:		
"A"	document defining the general state of the art which is not considered to be of particular relevance	"T" later document published after the international filing date or priority date and not in conflict with the application but cited to understand the principle or theory underlying the invention
"E"	earlier application or patent published on or after the international filing date	"X" document of particular relevance; the claimed invention cannot be considered novel or cannot be considered to involve an inventive step when the document is taken alone
"L"	document which may throw doubts on priority claim(s) or which is cited to establish the publication date of another citation or other special reason (as specified)	"Y" document of particular relevance; the claimed invention cannot be considered to involve an inventive step when the document is combined with one or more other such documents, such combination being obvious to a person skilled in the art
"O"	document referring to an oral disclosure, use, exhibition or other means	"&" document member of the same patent family
"P"	document published prior to the international filing date but later than the priority date claimed	
Date of the actual completion of the international search 25 November 2003 (25.11.2003)		Date of mailing of the international search report <b>18 DEC 2003</b>
Name and mailing address of the ISA/US Mail Stop PCT, Attn: ISA/US Commissioner for Patents P.O. Box 1450 Alexandria, Virginia 22313-1450 Facsimile No. (703)305-3230		Authorized officer <i>Felicia D. Roberts for</i> Tammy K. Field Telephone No. (703) 308-1235

# INTERNATIONAL SEARCH REPORT

PCT/US03/06715

## C. (Continuation) DOCUMENTS CONSIDERED TO BE RELEVANT

Category *	Citation of document, with indication, where appropriate, of the relevant passages	Relevant to claim No.
X	BOUVET et al. Evaluation of two colored latex kits, the Wellcolex Colour Salmonella test and the Wellcolex Colour Shigella test, for serological grouping of Salmonella and Shigella species. J. Clin. Microbiol. August 1992, Vol. 30, No. 8, pages 2184-2186, entire article.	1-2, 4, 9-10, 16, 18
---		-----
Y		3, 5-8, 11-15, 17